

TASKA CX™



# Sales Quick Reference Guide

## Introducing TASKA

TASKA hands instil confidence through aesthetics and reliability.

All hands are reliable, waterproof, easy to use.

Gen1/2 established our reputation as the #1 choice for specialists,  
CX enhanced it.

### Refined

1

Natural proportions and anatomical curves mimic the natural hand form.

2

The lightest multi-articulating hand with a flexion wrist.

3

Smooth, fast, responsive motions.

### Reliable

1

Patented shock-absorbing knuckle block, reduces impact-damage.

2

Knuckle-release system reduces damage from accidental overload.

3

IP67 waterproof without a glove.

### Empowering

1

Finger splay and Flexitool securely grips cutlery, tools and more.

2

Button panel and MyoTrigger control options make grip changing and locking intuitive.

3

Capacitive touch on all digits, including the thumb.

do more™

 **TASKA**™  
by eqwal.